

NEW

ASTM Colour Solutions

Manufactured by Reagecon

Reagecon
Delivering the Correct Result

Why use Reagecon's ASTM, Saybolt and Platinum-Cobalt Colour Solutions?

Commercial Benefits	Technical Benefits
<ul style="list-style-type: none">• Ready to use for use as calibration and/or quality control standards• Designed specifically for use in ASTM analytical methods• Presented in high quality tamper evident bottles• Customised standards available	<ul style="list-style-type: none">• Produced in accordance with ASTM methods (D1500, D6045, D1209)• Consistency of product – independent, traceable, certified.• Certificates of Analysis and SDS available online



Reagecon manufactures the full range of ASTM, Saybolt and Platinum-Cobalt Colour Solutions for use with ASTM standard methods including D1500, D6045, and D1209.

The products range from –

ASTM Colour Standard Sample A1- A7
Saybolt Colour Standards S+30 to S-15
Platinum-Cobalt Scale No. 0 – No. 500

These products are prepared gravimetrically on a weight/weight basis. Both solute and solvent are weighed on a balance calibrated by Reagecon engineers using OIML traceable weights. Reagecon holds ISO 17025 accreditation for weighing, (ref. 265C) and the resulting certificate is issued in accordance with the requirements of ISO/IEC 17025. The concentration of each standard is verified using a high performance top of the range calibrated spectrophotometer. The calibration of the spectrophotometer involves the use of high purity ISO Guide 34 accredited spectrophotometric standards.

Telephone
061 472622

Email
sales@reagecon.ie

Online sales
www.reagecon.com

NEW

ASTM Colour Solutions

Manufactured by Reagecon

Reagecon
Delivering the Correct Result

ASTM Colour Standards

Code	Description	ASTM Colour	Pack Size (ml)	US EPA Methods
ASTMA101	ASTM Colour Standard Sample A1	1.0	100	D6045, D1500
ASTMA105			500	
ASTMA301	ASTM Colour Standard Sample A3	3.0	100	
ASTMA305			500	
ASTMA501	ASTM Colour Standard Sample A5	5.0	100	
ASTMA505			500	
ASTMA701	ASTM Colour Standard Sample A7	7.0	100	
ASTMA705			500	

Saybolt Colour Standards

Code	Description	Saybolt Colour	Pack Size (ml)	US EPA Methods
SAYP301	Saybolt Colour +30	S+30	100	D6045
SAYP305			500	
SAYP251	Saybolt Colour +253	S+25	100	
SAYP255			500	
SAYP191	Saybolt Colour +19	S+19	100	
SAYP195			500	

Telephone
061 472622Email
sales@reagecon.ieOnline sales
www.reagecon.com

NEW

ASTM Colour Solutions

Manufactured by Reagecon

Reagecon
Delivering the Correct Result

Saybolt Colour Standards

(Continued)

Code	Description	Saybolt Colour	Pack Size (ml)	US EPA Methods
SAYP151	Saybolt Colour +15	S+15	100	
SAYP155			500	
SAYP121	Saybolt Colour +12	S+12	100	
SAYP125			500	
SAYP01	Saybolt Colour +0	S0	100	
SAYP05			500	
SAYN151	Saybolt Colour -15	S-15	100	
SAYN155			500	

Platinum-Cobalt Colour Standards *

Code	Description	Saybolt Colour	Pack Size (ml)	US EPA Methods
HAZ0	Platinum-Cobalt Colour 0	0	1000	D1209
HAZ10	Platinum-Cobalt Colour 10	10	1000	
HAZ25	Platinum-Cobalt Colour 25	25	1000	
HAZ40	Platinum-Cobalt Colour 40	40	1000	
HAZ50	Platinum-Cobalt Colour 50	50	1000	
HAZ80	Platinum-Cobalt Colour 80	80	1000	
HAZ100	Platinum-Cobalt Colour 100	100	1000	
HAZ250	Platinum-Cobalt Colour 250	250	1000	
HAZ500	500	500	1000	

* Standards with intermediate Platinum-Cobalt values are available on request

Telephone
061 472622Email
sales@reagecon.ieOnline sales
www.reagecon.com